

nulatex[®]

PROBE COVER

RELIABILITY REDEFINED



nulatex[®]

Nulatex Sdn Bhd is a privately owned manufacturing company of medical latex products based in Kluang, Johor, Malaysia.

Established in 2004, Nulatex Sdn Bhd produces and markets quality latex products to over 60 countries around the world. Located close to the suppliers of our key raw material, i.e. natural rubber latex, we are able to offer our customers superior products at competitive prices.

Lot 1870, 4th Mile, Jalan Mengkibol,
86000 Kluang, Johor, Malaysia.

+607 777 1801

info@nulatex.com

+607 777 1805

www.nulatex.com

@feelbynulatex

@feelbynulatex

Nulatex Sdn Bhd

Nulatex probe covers are manufactured from the highest quality natural rubber latex. Manufactured at an ISO 13485, CE, and US FDA listed manufacturing facility. Suitable for use on probes of Ultrasound Imaging System.



NULATEX PROBE COVER G (STERILE)

Individually sealed in square sachets.
 Single use (disposable), non-lubricated, Gamma Ray sterilized
 Diameter : 28mm / 34mm
 Length : 200 ± 5mm
 Width : 40 ± 2mm / 53mm ± 2mm
 Thickness : 0.056 - 0.065mm / 0.066 - 0.075mm



100 pcs



NULATEX PROBE COVER (NON STERILE)

Individually sealed in square sachets.
 Single use (disposable), non-lubricated (option of lubricated with silicone oil)
 MMDA Registration No: GA2997320-43686
 Diameter : 28mm / 34mm
 Length : 200 ± 5mm
 Width : 40 ± 2mm / 53mm ± 2mm
 Thickness : 0.056 - 0.065mm / 0.066 - 0.075mm



100 pcs



NULATEX PROBE COVER BULK PACK (NON STERILE)

Bulk naked, tested and rolled.
 Single use (disposable), non-lubricated
 MMDA Registration No: GA2997320-43686
 Diameter : 28mm / 34mm
 Length : 200 ± 5mm
 Width : 40 ± 2mm / 53mm ± 2mm
 Thickness : 0.056 - 0.065mm / 0.066 - 0.075mm



100 pcs